


VG 55 Ventilator

High-Flow Therapy (HFT)

- An air/oxygen blender, allowing for 21% to 100% FiO₂, generates up to 60L/min flow. HFT has been approved as an effective in the treatment of acute respiratory failure. Since 2011, the application of HFT has been increasing high in clinical field.



Complete Ventilation Modes for Clinicians

- Volume Guaranteed Pressure Support Ventilation (VGPS) assures tidal volume with auto adjusting pressure support, decreasing the workload for clinicians while ensuring patient safety.
- CPAP with Crelief option offers two levels of flow-based expiratory pressure relief, which improves patient comfort and treatment acceptance.



VG 55 Ventilator

CE 0123

AEO/MED
Reliable Quality Thoughtful Service

HQ: No.4 Hangfeng Road, Fengtai Science Park, Fengtai District, Beijing, China (100070)
Science Park: No. 10, Chaobai Street, Yanjiao Development Zone, Sanhe City,
Hebei Province, China (065201)
TEL: +86-10-5841 1198 FAX: +86-10-6371 8989
Http://www.aeonmed.com E-mail: marketing@aeonmed.com

Information contained in this leaflet is correct at the date of publication. The policy of Beijing Aeonmed Co., Ltd. is one of continued improvements to its products. Because of this policy Beijing Aeonmed Co., Ltd. reserves the right to make any changes, which may affect the information in this leaflet without giving prior notices.
AVT60-1801

AEO/MED
Reliable Quality Thoughtful Service

VG 55 Ventilator

Start Your Natural Breath with VG55

VG55 is specially designed for Non-Invasive Ventilation (NIV). Research has shown that early application of NIV can effectively decrease intubation rates and hospitalization length. The application of NIV can also prevent ventilator associated pneumonia (VAP), as well as promote respiratory care. VG55 adopts advanced turbine technology; provides NIV with ICU-grade monitoring.



Standby mode for patient-clinician interaction

Various information interfaces:
USB, VGA, Nurse Call,
Wifi & Bluetooth

Smart trigger and cycling ensures patient-machine synchrony

Alarm reset helps maximize patient safety while reducing excessive alarms

15 Seconds

Start Your Natural Breath within 15 Seconds

- Auto-calibration is performed within 15 seconds, safely guaranteeing the ventilator's performance.
- Quick setup with guide-way system and intuitive interface.
- Smart trigger combines pressure and flow triggers, highly improving the patient-machine synchrony, allowing the patient to breathe in a more natural manner.

5 Minutes

Become a Respiratory Therapist in 5 Minutes

- With 15.6" intelligent touchscreen, the VG55 displays complete ventilation monitoring in a single screen.
- Simplifies the patient management with user-friendly interface and guide-way information.
- Brief explanations of all included ventilation modes provide a clear clinical reference for treatment options.

in Seconds

Get What You Need in Seconds

- Two minutes of emergency oxygen therapy provides 100% oxygen for whenever and wherever needed.
- Extended SpO₂, PR and EtCO₂ module fulfill clinical requirements.
- Advanced integrated Aerogen nebulizer enhances respiratory care.
- High-Flow Therapy up to 60L/min.

